Schmerzlinderung Pain Relief

durch elektrische Nervenstimulation (TENS)
through electrical nerves stimulation (TENS)

Lieferant / distributor:

Hersteller / manufacturer:

Endelner Feld 9 • D - 46286 Dorsten
www.saneo-line.com
info@saneo-line.com
Hotline +49 (0) 18 05 . 18 15 18
(0,14 €/Min. aus dem deutschen Festnetz,
Mobilfunk max. 0,42 €/Min.) (0,14 €/min. from
a German landline, other providers may differ)

05/2010 / R1 CG

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Notizen:

Dear customer ...

Many thanks for trusting in us and selecting this top-quality product from

tic Medizintechnik GmbH & Co. KG.

These instructions for use will guide you through the application options of the therapy device. In addition, we provide you with therapy tips.

Please always keep the instructions for use with the device.

en

This device does not replace a visit to the doctor!

Although application of the therapy device may achieve successful results, a significant reduction in pain should not be equated with cure of the root of the pain.

Device overview

- 19 therapy programs
- 3 configurable programs
- Timer function
- Large LC display: full and clear representation of all functions and parameters
- Automatic switch-off after treatment
- Simple operation
- Battery monitor
- Automatic electrode recognition
- No-load function
- Belt clip for greater mobility

Contents

Important notes	3
... on these instructions for use	3
... on the specific function	3
... on liability and warranty	3
Safety notes	4
Device description	6
TENS pain therapy using the SaneoTENS	6
Contra-indications / application limitations	7
Device diagrams	8
Scope of delivery	10
Accessories	10
Operating the SaneoTENS	11
Possible displays	11
Preparing for treatment	12
Positioning the electrodes	12
Switching the SaneoTENS on and off	13
Setting the program	13
Selecting the stimulation intensity	14
Ending treatment & cleaning the therapy system	14
Exchanging the battery	15
Program overview	16
... according to program	16
... according to indications / body regions	18
Programming unrestricted-use programs	20
Parameter overview	20
Entering parameters	21
Technical data and support	22
Symbols	22
Error and problem handling	23

Important notes

en

... on these instructions for use

Please read these instructions for use prior to initial use.

These instructions for use are aimed at doctors and therapists as well as patients. They form part of the device and must be available to the user/patient at all times. Excerpts are not permitted.

Please carefully read the notes on proper use (page 3), the Safety notes (page 4) and the contra-indications (page 7).

... on the specific function

SaneoTENS is a home therapy device, which serves to provide electrical nerve stimulation through the skin. Please use it exclusively to treat indications described in the chapter Program overview (page 18).

Any application not described in the instructions for use, or described as impermissible, shall be considered to be improper use of the device.

Only accessories authorised by the manufacturer may be used. (e.g. TENS body surface electrodes).

... on liability and warranty

The warranty applies 24 Months to the SaneoTENS therapy device and all accessories (excluding consumable supplies) from the day of delivery. During this period of time the device/accessories can be repaired or replaced by the manufacturer at no charge.

The manufacturer assumes no liability for damages resulting from improper use, non-adherence to these instructions for use and unauthorised modifications to the system. The warranty shall become null and void.

The manufacturer has the right to perform technical changes, which approve the therapy device, his functions and the accessories.

Important notes

Safety notes

SaneoTENS

is a medical device of class IIa (MDD). When used commercially, a safety-related check must be performed by the manufacturer every two years.

SaneoTENS

- Use only with original accessories
- Protect from water and other liquid substances and do not use at humidity exceeding 90%
- Use only at temperatures between 10°C and 40°C
- Store only at temperatures between 0°C and 50°C
- Do not use if device is faulty or damaged (e.g. due to a fall)
- Exercise particular care and attention when using the device on children
- Ensure that children and unauthorised persons are unable to access the device
- Store device in packaging and protect from damage and contamination
- Do not use while driving or operating machinery
- Do not use when battery case is open

Important notes

CAUTION

Interaction with other devices/systems

Interaction with other (medical) devices may impact the output values. Do not use the SaneoTENS near:

- Microwave or short-wave devices
 - Portable HF communication devices (mobile phones)
- ➔ During treatment, please ensure that your mobile phone is at least three metres away.

en

CAUTION

System fault due to break in electrode cable

Improper use, extreme bending or irregular pulling of a damaged electrode cable can lead to system faults.

- ➔ Do not force the electrode cable into a bent position.
- ➔ Do not pull the cables.

CAUTION

Electro-chemical burns!

Do not use SaneoTens when simultaneously connecting to HF surgical or electro-coagulation devices. Disregard may result in electro-chemical burns.

CAUTION

Uncontrollable current flow!

Touching metallic parts inside of the battery case can cause uncontrollable discharges. Use only with closed battery case!

Device description

TENS pain therapy using the SaneoTENS

TENS (Transcutaneous Electrical Nerve Stimulation) has been medically recognised and proven for over 40 years.

SaneoTENS is a TENS therapy system for the effective treatment of acute and chronic pain.

During TENS therapy a medically effective current that is non-hazardous to health is directed through the skin to the tissue via two or four electrodes.

Two frequency ranges are available in the treatment programs.

- TENS therapy (1 Hz to 20 Hz)

... is ideally suited for treatment of chronic pain as it can stimulate the release of pain-reducing substances that are produced naturally in the body (e.g. endorphins, dynorphins). The body's own painkillers also remain effective after treatment.

- TENS therapy (50 Hz to 100 Hz)

... is particularly suited for treatment of acute pain. It can block or reduce pain signals to the brain.

Treatment frequency

Where necessary, you can perform several treatments daily, provided no side effects such as skin reactions occur.

Device description

en

Contra-indications / application limitations

Medical circumstances under which the SaneoTENS may not be used:


- During pregnancy, the device may never be used in the abdominal area. Current used for medical purposes can lead to miscarriage!
- When operating machinery (e.g. driving a car) or while sleeping
- In case of metallic implants in the area of treatment
- In the event of skin defects or defects of the mucous membrane (injuries) and acute infections of the mucous membrane, skin or subcutaneous skin
- In case of cardiac arrhythmia or heart muscle disease
- Where psychological limitations and/or cognitive differences are present
- Never use in the heart region or near the carotid artery
- Caution! The application on the thorax may lead to ventricular fibrillation.
- The device is not suitable for persons with a pacemaker and/or active implants (defibrillators)! 
- Current used for medical purposes promotes blood circulation. Treatment must be stopped and/or performed under medical supervision if:
 - A drastic change is experienced in the sensation of the current, even though the stimulation intensity was not changed.
 - Previously unknown infections or bleeding is noticed in the area of treatment prior to commencing treatment.
 - Increasing or strong skin changes are noticed, in which case a doctor must be consulted (possible allergy).

Device description

Device diagrams

Overall view

en


- 1** Electrode connector socket, Channel A
- 2** Electrode connector socket, Channel B
- 3** ▲: Increase intensity
- 4** ▼: Decrease intensity
- 5** **Menu** button (select program)
- 6** **On/Off** button

Device description

Displays

All displays are listed in the chapter Possible displays, page 11.

en

- | | | | |
|----------|---|-----------|---------------------------------|
| 1 | Battery charge status | 7 | Stimulation intensity Channel B |
| 2 | Frequency modulation active | 8 | Stimulation intensity Channel A |
| 3 | Current stimulation phase (flashing) | 9 | Burst mode active |
| 4 | Remaining treatment time in minutes:seconds | 10 | Invalid entry by active locking |
| 5 | Locking active | 11 | Set program |
| 6 | Therapy session has been interrupted. | | |

Device description

Scope of delivery

Item no.	Designation
ti2003	1 x SaneoTENS therapy device 1 x Instructions for use
ti6060	1 x Block battery (9 V - alkali)
ti2013	1 x TENS body surface electrodes 50 mm x 50 mm (4 units)
ti2033	2 x Cable for connections to body surface electrodes (2-pole), 2 mm

Accessories

Item no.	Designation
	TENS body surface electrodes:
ti2011	32 mm round (4 units), PZN 3888368
ti2013	50 x 50 mm (4 units), PZN 3888546
ti2015	50 x 90 mm (4 units), PZN 3888351

In Germany, selected items are also available from a pharmacy. In this case, please provide the PZN (pharmaceutical central number).

Accessory diagrams

TENS body surface electrodes

32 mm round
(ti2011)

50 mm x 50 mm
(ti2013)

50 mm x 90 mm
(ti2015)

Besides the items listed above, you can also find additional SaneoLINE products at our website:

<http://www.saneo-line.com>

Operating the SaneoTENS

en

Possible displays

Pxx

Treatment program (active program)

Locking. The keypad is locked five seconds after the last entry to prevent unintentional entries. Unlocking only by means of keystroke ▼.

Invalid entry / invalid keystroke

Locking remains active.

Therapy session interrupted due to insufficient patient contact.

See chapter Error and problem handling (page 23), item 2.

Treatment time per session in minutes (**min**) and seconds.

mA

Stimulation intensity of 0 mA to 60 mA.

Burst

Mode for preventing acquired tolerance tendencies.

FM

(**FM**) Frequency modulation

Units of displayed values:

Hz

Frequency (pulses per minute)

µs

Pulse duration (in microseconds)

sec

Stimulation and/or resting phase (in seconds)

min

Treatment duration (in minutes)

Battery charge status:

The battery is charged; the device is ready for use.

The battery and the device are ready for use.

Flashing: Battery must be replaced in the short term; battery charged state critical.

The current stimulation and/or resting phase is indicated by flashing. Image from left:


Pause | ▲Increase | ■Plateau | ▼Descent | Pause

Operating the SaneoTENS

Preparing for treatment

The device is switched off; no display is present.

- Please connect the electrode cable with connector socket Channel A at the front of the device.
- Should you intend using both channels, connect an additional electrode cable with connector socket Channel B.

Preparing skin areas to be treated

Clean the skin areas to be treated of all perspiration, moisture and dirt. Ensure that the skin is dry (do not use any ointments or creams) and intact.

- Please remove the body surface electrodes from the protective packaging.
- Connect the electrodes with the electrode cable (scope of delivery).

Positioning the electrodes

Electrode application can be used to treat various indications. An electrode with two adhesive surfaces is connected to Channel A or to Channel B.

The electrode application and the program selected must match; see chapter Program overview, page 18.

- Place (stick) the electrodes as appropriate to match your program and according to the electrode application illustrated **in the cover of these instructions for use**. Observe the safety notes (page 4)!
- After treatment, the electrodes can be stuck back on the carrier sheet and stored in the protective packaging.

Operating the SaneoTENS

Switching the SaneoTENS on and off

Use the **On/Off** button to switch the device on and off. The software version is briefly displayed after switch-on.

en

Setting the program

Set your required program as follows (see chapter Program overview, page 16):

- ➔ Please switch the therapy device on by means of the **On/Off** button. The *treatment menu* of the preset program then briefly appears in the display.
- ➔ Hold the **Menu** button down.
After approx. three seconds the *program selection menu* is displayed with the current program.
- ➔ Press the **▲** button several times until the required program is displayed.
- ➔ Briefly press the **Menu** button.
The *program selection menu* closes. The selected program is immediately displayed in the *treatment menu*.

Now please select the individual stimulation intensity.

Operating the SaneoTENS

Selecting the stimulation intensity

en

Take note of your own perceptions!

You should clearly feel the current or stimulation intensity, but it should not be an unpleasant sensation and certainly not painful.

The stimulation intensity for Channel A and/or Channel B is displayed on the left and/or right-hand side of the *treatment menu* (unit: **mA**).

Please select the stimulation intensity for **Channel A** using the ▲/▼ buttons to the **left** of the display and for **Channel B** using the ▲/▼ buttons to the **right** of the display.

- ➔ Pressing the ▲ button briefly increases the stimulation intensity by 1 mA for the relevant channel.
- ➔ Pressing the ▼ button briefly decreases the stimulation intensity by 1 mA for the relevant channel.

Treatment is started.

Where the appropriate stimulation intensity was set and not changed for a period of 15 seconds, the **Lock** symbol  appears in the display. The keypad is locked.

- ➔ Where required, unlock the keypad by means of the ▼ button, e.g. to adjust the intensity.

Ending treatment & cleaning the therapy system

The device automatically switches off when the treatment time has ended. Treatment can be aborted at any time by using the **On/Off** button.

- ➔ Please do not remove the body surface electrodes before the device has been switched off.
- ➔ Use a damp, lint-free cloth to clean the electrodes after use. This will increase the durability of the electrodes.

Operating the SaneoTENS


- ➔ Stick the electrodes on the carrier sheet and store the carrier sheet in the protective packaging.

If the adhesiveness of the electrodes decreases, please use water to wet the adhesive surfaces. If this is no longer effective, please replace the electrodes.

en

Exchanging the battery

- ➔ Please switch the device off before exchanging the battery.

- 1** Open the battery compartment.
 - 2** Place the battery in the compartment in the correct position.
 - 3** Close the battery compartment.
- ➔ Please take note of the national statutory regulations applicable in your country in terms of environmentally friendly disposal of batteries, accessories and packaging.

Program overview

... according to program

The following table lists all programs (**Prog**), their **application area**, **mode of operation**, **settings** and **duration**.

Prog	Application area	Mode of operation	Frequency (Hz)	Pulse width (μ s)	Rise time (sec)	Plateau time (sec)	Descent time (sec)	Pause (sec)	Duration (min)
1	Acute, neuropathic or chronic pain	continuous stimulation, Channel A and B identical	100	200					30
2	Acute, neuropathic or chronic pain	cyclic stimulation, Channel A and B identical, 3 sec 100 Hz, 3 sec 3 Hz	100/3	180					30
3	Alternative to Program 1 and 2, especially for chronic pain	continuous stimulation, Channel A Channel B (attach Channel A and B crosswise)	100 4	200 200					30
4	Alternative to Program 1, 2 and 3 to prevent therapy resistance	Burst stimulation, Channel A and B identical, 7 pulses per burst group, pulse frequency 2 Hz, burst frequency 100 Hz	100	250					30
5	Alternative to Program 1, 2 and 3 to prevent therapy resistance	frequency modulation, Channel A and B identical	8-100	250	2	2	2	2	30
6	Muscle pain	continuous stimulation, Channel A and B identical, 10 min 100 Hz, then 10 min 4 Hz then 10 min 100 Hz,	100/4	150@ 100Hz 200@ 4HZ					30
7	Muscle tension	cyclic stimulation, Channel A and B identical, 30 sec 4 Hz, 30 sec 6 Hz, 30 sec 8 Hz, 20 sec 4 Hz, 20 sec 6 Hz, 20 sec 8 Hz in between the phases 2 sec pause	4/6/8	250					30
8	Muscle tension	amplitude modulation, Channel A and B alternating	75	150	1	0	1	2	30

Program overview

en

Prog	Application area	Mode of operation	Frequency (Hz)	Pulse width (μ s)	Rise time (sec)	Plateau time (sec)	Descent time (sec)	Pause (sec)	Duration (min)
9	Muscle tension	cyclic stimulation, Channel A and B identical, 5 sec 40 Hz, 10 sec 4 Hz, 5 sec 50 Hz	40/4/5 0	250					30
10	Spine syndromes (cervical spine, vertebral, lumbar spine)	continuous stimulation, Channel A and B identical	110	180					30
11	Spine syndromes (cervical spine, vertebral, lumbar spine)	continuous stimulation, Channel A and B identical	80	150					30
12	Knee, joint or tendon complaints, phantom pain	continuous stimulation, Channel A and B identical	40	200					30
13	Rheumatic joint complaints	frequency modulation, Channel A and B identical	2-8	250	2	2	2	2	30
14	Neuropathic pains, shingles	continuous stimulation, Channel A Channel B	10 2	200 200					30
15	Menstrual complaints	continuous stimulation, Channel A and B identical	8	200					30
16	Blood circulation disorders	continuous stimulation, Channel A and B identical	2	250					30
U1	Free, see chapter Programming unrestricted-use programs	basic program continuous TENS stimulation, frequency and duration configurable	2-120	250					5-99
U2	Free, see chapter Programming unrestricted-use programs	basic program frequency modulation, frequency, rise, plateau, descent and duration configurable	2-100	250	0,5-5,0	1-10	0,5-5,0	1-10	5-99
U3	Free, see chapter Programming unrestricted-use programs	basic program burst stimulation, Channel A and B identical, 7 pulses per burst group, pulse frequency 2 Hz, burst group intervals and duration configurable	100	250				0,1-5,0	5-99

Program overview

... according to indications / body regions

The programs 1 to 5 basically can be used for treating almost every kind of pain. Furthermore the following table lists alphabetically the most common indications resp. the body regions where pain appears most commonly. It gives you recommendations for accomplishing therapy with further special programs and for electrode applications (see the cover of this instruction manual).

To receive maximum results from your stimulation, it is important to place the electrodes correctly. The optimum placement varies slightly from person to person, so try moving the electrodes around until you get the most comfortable and effective results. Before placing the electrodes on your body, be sure to always wash and dry the skin in order to ensure good conductivity. Be sure to place electrodes only on healthy skin. Furthermore you always should place the electrodes without overlapping.

Indication / body region	Program	Electrode application
Achilles tendon	12	40
Ankle	12	38
Blood circulation disorders	16	28
Buttocks musculature	6	29
Elbow	12	24
Forearm (tennis elbow)	6	25, 26, 27
Hip	12	30, 31
Knee joint	12	35
Lower leg	6	36, 37
Shingles	14	*

Program overview

en

Indication / body region	Program	Electrode application
Sciatic nerve	14	15, 16, 17
Menstrual complaints	15	34
Migraine	1	3
Neck and back	6	2, 4, 5, 6, 7, 8, 9, 11
Rheumatic joint complaints	13	*
Shoulder	6, 12	10, 12, 39
Spine (Cervical spine, vertebral, lumbar spine syndrome)	10, 11	5, 11, 14, 15, 16, 18
Stomach musculature	6	19, 20, 21
Tension headaches	7, 8	2
Tensions	7, 8, 9	*
Thigh	6	32, 33
Trigeminal neuralgia	14	1a, 1b, 1c
Upper arm	6	22, 23

* Place the electrodes on or around the pain area.

Program overview

Programming unrestricted-use programs

Therapy programs **U1**, **U2** and **U3** can be changed as required.

U1: Basic program – continual stimulation

U2: Basic program – frequency modulation

U3: Basic program – burst stimulation

Parameter overview

Depending on the basic program, only specific parameters can be set. Please refer to the table below for sequence and significance of the individual parameters.

The value range of the parameters is indicated (in brackets).

Programming of medically effective and meaningful parameter values that can be combined should only occur after consultation with your doctor or therapist.

Parameter	Description (value range)
Hz	Treatment frequency in Hz (2 to 120)
	Rise time in sec (0.5 to 5)
	Plateau time in sec (1 to 10)
	Descent time in sec (0.5 to 5)
	Pause time in sec (1 to 10)
dEL	Pause between burst groups in sec (0.1 to 5)
	Treatment time in min (5 to 99)

Program overview

en

Entering parameters

You can enter settings for therapy programs **U1**, **U2** or **U3**. Your program configuration will then be called up as for a default program.

- ➔ Hold the **Menu** button down.
After approx. three seconds the *program selection menu* is displayed with the current program. 
- ➔ Press the ▲ button several times until the therapy program to be configured is displayed.
- ➔ Press the **Menu** button to navigate to the first user-defined parameter.
The parameter value flashes and can be reduced by actuating the ▼ button or increased by actuating the ▲ button.
- ➔ Please press the **Menu** button to confirm the value selected.
The next parameter is displayed.
- ➔ Please continue as described until you reach the last parameter (**treatment time**).

Once you have also confirmed the treatment time using the **Menu** button, the *program selection menu* closes.

The *treatment menu* for the configured program appears in the display. Treatment with the preset parameter values can now begin.

Starting the treatment

To start the treatment, please set the stimulation intensity as described in the chapter *Selecting the stimulation intensity*, page 14.

Technical data and support

Technical data

en

Name and model	SaneoTENS Pain Relief
Type	2-channel electro-stimulation system
Dimensions	105 mm x 65 mm x 26 mm
Weight	Including battery approx. 130 g
Voltage supply	Alkaline block battery (9 V) 6LR61
Classification	BF 
Current pulse shape	Biphasic rectangular pulse, charge-compensated
Relative humidity	10% to 90%
Operating temperature	0°C to 40°C
Storage temperature	0°C to 50°C
Channels	2
Current intensity	Constant Current, max. 60 mA at 500 Ω (± max. 5 mA)
Maximum output voltage	100 V ±15%
Treatment time	Can be adjusted from 5 min. to 99 min.
Pulse width	200 µs to 250 µs
Frequency range	2 Hz to 120 Hz

Symbols

Floating application component, **Body Floating** type. The product provides protection against electrical shock in the case of permissible discharge current.

Proof of conformity in terms of Directive 93/42/EEC of the Board dated 14 June 1993 regarding medical devices. This is a medical device of class IIa (MDD).

Year of manufacture

This product carries the recycling symbol in accordance with EU Directive 2002/96/EC. Disposal via domestic or residual waste is not permitted. The dealer / manufacturer is the contact person responsible for disposal.

Device serial number

Lot designation

Usable until

Error and problem handling

In the event that the recommended solutions do not solve the problem, please contact Customer Services or your dealer.

1. The device does not switch on when the On/Off button is actuated.

- Ensure that the battery has been inserted correctly and contact has been established.
- Replace battery, where necessary.

2. No stimulation can be felt.

- Press the **On/Off** button to start the device.
- Check the position of the electrodes. The electrodes may not overlap.
- Gradually increase the stimulation intensity.
- Check the cable connections between the device and the electrodes.
- Please restore the adhesiveness of the electrodes by moistening.

If the adhesive surfaces of the electrodes are depleted, current distribution is no longer guaranteed.

Please replace the electrodes.

Notes: