Muskeltraining Muscle Training

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Dear customer ...

Many thanks for trusting in us and selecting this top-quality product from

tic Medizintechnik GmbH & Co. KG.

These instructions for use will guide you through the application options of the muscle training device. In addition, we provide you with application tips.

Please always keep the instructions for use with the device.

en

Device overview

- 20 training programs
- 2 configurable programs
- Timer function
- Large LC display: full and clear representation of all functions and parameters
- Automatic switch-off after training
- Simple operation
- Battery monitor
- Automatic electrode recognition
- No-load function
- Belt clip for greater mobility

Contents

Important notes	3
... on these instructions for use	3
... on the specific function	3
... on liability and warranty	4
Safety notes	5
Device description	7
The SaneoSPORT muscle training device	7
Program summary	8
Preparation for training session	9
1. Warm up I	9
2. Warm up II	9
Muscle build-up training	9
Support for active sports training	10
9. Aerobic endurance training	10
10. Anaerobic endurance training	10
11. Strength and speed training	11
12. Explosive strength training	11
13. Crosstraining (overall conditioning)	12
Training follow-up	12
14. Relaxation	12
15. New strength	12
16. Stretching	12
Pain relief	13
Configurable programs	13
Contra-indications / application limitations	14
Scope of delivery	15
Accessories	15
Device diagrams	16
Possible displays	18
Operating the SaneoSPORT	19
Preparing for application	19
Positioning the electrodes	19
Switching the SaneoSPORT on and off	20
Setting the program	20
Selecting the stimulation intensity	21
Programming configurable programs	22
Parameter overview	22
Entering parameters	23
Ending training & cleaning the electrodes	23
Exchanging the battery	24
Technical data and support	25
Symbols	25
Error and problem handling	26

Important notes

en

... on these instructions for use

Please read these instructions for use prior to initial use.

These instructions for use form part of the device and must be available to the user at all times. Excerpts are not permitted.

Please carefully read the notes on proper use (page 3), the Safety notes (page 5) and the Contra-indications (page 14).

... on the specific function

The SaneoSPORT is a home treatment device for the stimulation of various kinds of nerve fibres.

1. Motor nerves, in order to stimulate muscles to contract or relax. This is generally described as electrical muscle stimulation. Benefits depend on the stimulation parameters.

According to program parameters (impulses per second, pulse duration, pause time etc.) the SaneoSPORT serves for warming up muscles, muscle rehabilitation, blood circulation promotion or for solving muscle tensions. By appealing the different muscle fiber kinds (slow, intermediate and quick fibers) by means of various contraction speeds the SaneoSPORT helps to train the different kinds of muscle work. For example a cross-country skier has predominantly slow fibers while with the sprinter the quick fibers prevail.

2. Certain sensitive nerve fibres, in order to achieve a pain-relieving effect.

According to kind of current the pain transference to the brain is blocked or the sensation of pain is decreased by an increased endorphin production.

Important notes

Electro-stimulation with the SaneoSPORT can bring the following benefits for the user.

- The positive effect exerted on muscle characteristics is significant.
- Very little strain is placed on the joints and tendons.
- Resultant fatigue of the cardiovascular system is absent.
- Muscles can be required to take on a higher workload.

The greatest amount of benefit is achieved when electro-stimulation with the SaneoSPORT is supplemented by:

- healthy and balanced nutrition,
- regular movement,
- a balanced lifestyle

Any application not described in the instructions for use, or described as impermissible, shall be considered to be improper use of the device.

Only accessories authorised by the manufacturer may be used. (e.g. body surface electrodes).

... on liability and warranty

The warranty applies 24 Months to the SaneoSPORT muscle training device and all accessories (excluding consumable supplies) from the day of delivery. During this period of time the device/accessories can be repaired or replaced by the manufacturer at no charge.

The manufacturer assumes no liability for damages resulting from improper use, non-adherence to these instructions for use and unauthorised modifications to the system. The warranty shall become null and void.

The manufacturer has the right to perform technical changes, which approve the muscle training device SaneoSPORT, its functions and the accessories.

Safety notes

SaneoSPORT

is a medical device of class IIa (MDD). When used commercially, a safety-related check must be performed by the manufacturer every two years.

SaneoSPORT

- Use only with original accessories
- Protect from water and other liquid substances and do not use at humidity exceeding 90%
- Use only at temperatures between 10°C and 40°C
- Store only at temperatures between 0°C and 50°C
- Do not use if device is faulty or damaged (e.g. due to a fall)
- Exercise particular care and attention when using the device on children
- Ensure that children and unauthorised persons are unable to access the device
- Store device in packaging and protect from damage and contamination
- Do not use while driving or operating machinery
- Do not use when battery case is open

Important notes

CAUTION

Interaction with other devices/systems

Interaction with other (medical) devices may impact the output values. Do not use the SaneoSPORT near:

- Microwave or short-wave devices
 - Portable HF communication devices (mobile phones)
- ➔ During treatment, please ensure that your mobile phone is at least three metres away.

CAUTION

System fault due to break in electrode cable

Improper use, extreme bending or irregular pulling of a damaged electrode cable can lead to system faults.

- ➔ Do not force the electrode cable into a bent position.
- ➔ Do not pull the cables.

CAUTION

Uncontrollable current flow!

Touching metallic parts inside of the battery case can cause uncontrollable discharges. Use only with closed battery case!

Device description

en

The SaneoSPORT Muscle Training device

The SaneoSPORT muscle training device works with two forms of therapy which have been recognised in medical science for over 40 years:

EMS (Electrical Muscle Stimulation) and

TENS (Transcutaneous Electrical Nerve Stimulation).

EMS, which is an extremely effective therapy which is predominantly applied in rehabilitation, is a method which uses electrical impulses to address certain muscle groups in a targeted way. Its areas of deployment include the avoidance of muscle wastage and the re-establishment and relaxing of muscles.

Muscle activity is normally triggered by a transmission command in the brain. A pulse in the form of an electrical signal is sent to the responsible nerve pathway via the spine. This pathway supplies the nerve fibres of the muscle, and the pulse thus ultimately causes the muscle to tension. EMS generates artificial stimulation of the muscle by means of low-frequency electrical pulses that stimulate the nerves located beneath the skin, which in turn activate the corresponding muscle groups. The muscle is unable to distinguish between an impulse from the brain or one that comes from outside. As a result, it reacts to the stimulation in its normal way, tensions and consequently burns up the energy spent for the activity.

TENS is used for the effective treatment of acute and chronic pain. It involves electrical impulses which, depending on the stimulation frequency, either bring about a distribution of the body's own pain-relieving substances or else block transmission of pain to the brain.

The current which is conducted through the skin into the body tissues by the electrodes during the massage is medically effective and represents no health hazard.

Application frequency

Where necessary, you can perform several applications daily, provided no side effects such as skin reactions occur.

Device description

Program summary

The table below gives a summary of the SaneoSPORT muscle training programs and provides you with recommendations on the placement of the electrodes. The following pages contain detailed information on the individual programs.

Prog	Type	Description	Elektrode placements
1	Preparation for training	Warm up I	see Prog 3 - 8
2		Warm up II	see Prog 3 - 8
3	Muscle build-up training	Chest / shoulder and upper back musculature	2 - 7
4		Front and back upper and lower arm musculature (incl. biceps, triceps)	8 - 13
5		Straight and oblique abs musculature, lower back musculature	14 - 15, 17 - 19
6		Buttocks musculature	20 - 22
7		Thigh musculature	23 - 25
8		Lower leg musculature	27 - 29
9	Support for active sports training	Aerobic endurance training	see Prog 3 - 8
10		Anaerobic endurance training	see Prog 3 - 8
11		Strength and speed training	see Prog 3 - 8
12		Explosive strength training	see Prog 3 - 8
13		Crosstraining	see Prog 3 - 8
14	Training follow-up	Relaxation	see Prog 3 - 8
15		New strength	see Prog 3 - 8
16		Stretching	see Prog 3 - 8
17	Pain Relief	Pain Relief I	1 - 31
18		Pain Relief II	1 - 31
U1	Configurable programs	Muscle build-up training I	see Prog 3 - 8
U2		Muscle build-up training II	see Prog 3 - 8

Preparation for training session

1. Warm up I

... optimally prepares the muscles for unfamiliar sporting activities. Activation of the blood circulation and stimulation of the metabolism significantly improves the contraction capability of the musculature. An invigorating effect can generally already be felt during the application period in the form of a pleasant heat sensation. The **Warm up I** program should be used 30 minutes prior to the start of training.

2. Warm up II

... causes a dense capillary network to be formed around the muscle fibres. Capillaries are tiny blood vessels that enable a metabolic exchange to take place between the blood and the tissue. The denser the capillary network, the larger the exchange area. In this manner, the **Warm up II** program causes extreme intensification of the blood circulation and oxygen supply to the stimulated muscles. The glucose supply is also improved and the elimination of lactic acid, which increases when the muscles are extremely stressed, is accelerated. After applying this program you will experience any physical exertion or sporting activity with an enhanced feeling of well-being and your recovery time after the training session will be considerably reduced. **Warm up II** is particularly suited for persons who are only marginally fit or not fit at all, as preparation for all endurance sports. Application should take place regularly over several weeks prior to the start of training as well as between the actual training sessions.

Muscle build-up training

... adapted to the specific muscle characteristics of the individual body parts

3. **Chest / shoulder musculature, upper back musculature**
4. **Front and back upper and lower arm musculature (incl. biceps, triceps)**
5. **Straight and oblique abs musculature, lower back musculature**

Device description

6. **Buttocks musculature**
7. **Thigh musculature**
8. **Lower leg musculature**

Support for active sports training

9. Aerobic endurance training

Aerobic sports types include hiking, walking, swimming, cycling, jogging, cross-country skiing, ice skating, inline skating.

Endurance training is healthy for the heart and cardiovascular system. Optimal endurance training occurs in the aerobic range, when sufficient oxygen is available to the body and the oxygen intake and oxygen consumption are in equilibrium. As a result, the body can burn up fat. In addition, it only makes use of reserves that are instantly available to it, such as sugar.

General aerobic endurance is understood to mean training at a rate in which a conversation is still possible while walking, for example. The better the aerobic endurance (also known as basic endurance or general endurance), the quicker the recovery rate after training, and the sooner you will be ready for your next bout of activity.

The **aerobic endurance training** program considerably increases oxygen consumption of the muscles in the area of stimulation, and consequently trains and increases your oxygen intake and absorptive capacity. In this, tendons and joints experience almost no stress at all.

Note: You can achieve a particularly effective result by regularly combining this program with the **Warm up II** program.

10. Anaerobic endurance training

Anaerobic sports types include weight training, sprinting or stop-and-go sports types such as tennis and squash. The objective of anaerobic endurance training is to be able to sustain intensive exertion for as long as possible. In this, the body is not supplied with sufficient oxygen for energy production, resulting in the metabolism running at an oxygen debt. This undersupply of oxygen means that the body has to generate energy by different ways and means, e.g. by burning carbohydrates.

Device description

en

Producing energy in this manner results in an increased formation of lactic acid, and as the body is able to break down only small amounts of lactic acid, muscles can become over-acidic and fatigued if too much lactic acid is produced.

The **Anaerobic endurance training** program allows you to improve your performance capacity in the above-mentioned sports types, for example. It ensures greater tolerance of increased lactic acid concentrations and improves muscle resistance to fatigue.

Note: We recommend using the **Relaxation** program after particularly intensive training.

11. Strength and speed training

The following sports types require both strength and speed: competitive cycling, time-based swimming, team sports such as football, handball and volleyball. This program allows you to maximise the strength as well as the contraction speed of the stimulated muscles through the unique type of muscle contraction and the extreme effort required.

Note: By subsequent use of the **Relaxation** program, you can enhance the effect achieved through training even more.

12. Explosive strength training

Explosive strength training is typical for the following sports types: sprinting and short-distance swimming, high jump and long jump, hurdling, throwing sports etc. These sports types are generally trained for the purpose of competing. The preparation phase is consequently crucial. The capability of a muscle to develop a high level of strength in the fastest time possible is known as muscular explosive strength. As a rule, tough active training is a prerequisite for good results in these sports types. More frequently than not, excessive load is placed on the affected muscles, which tire quickly as a result. Regular use of the **Explosive strength training** program will enable you to significantly improve your explosive strength performance while simultaneously reducing the number of active training sessions.

Device description

13. Crosstraining (overall conditioning)

The **Crosstraining** program comprises different work-out sequences and is aimed especially at recreational sportspeople who want to achieve a good level of fitness in all types of muscle activity (aerobic and anaerobic endurance, strength, speed and explosive strength). High-performance athletes should however also use the **Crosstraining** program at the start of the season or after longer periods of inactivity to get their muscles used to all the different types of load again. The training sessions can subsequently be specified and increased as required.

Training follow-up

14. Relaxation

The **Relaxation** program is recommended for use after training and competitive sport. It is based on medical applications and aims to release tense muscles, increase the blood circulation and relieve pain. This program helps regenerate your muscles more quickly after experiencing intensive loads. The program should be used within 3 hours of training and should be supported by a hot shower and large liquid intake, for example. You will already start feeling that pleasant relaxing sensation during application. You can intensify the effect by a follow-up treatment after approximately 6 to 8 hours, where muscles have been particularly worked.

15. New strength

The **New strength** program is perfect as a supplement or replacement for an active training session to restore your muscles "on the day after". Light aerobic muscle training in the form of light, non-tiring, continued contraction is accompanied by pain relief and increased blood circulation.

16. Stretching

The **Stretching** program works on the physiological principle of reciprocal innate reflexes and causes distinctive loosening of the musculature. Active stretch exercises should be supplemented by placing the electrodes on the opposing muscle. This results in an intensified stretch effect.

Pain relief

The programs **17.** and **18.** are optimal for treating chronic and acute sports-induced pain, e.g. sore muscles, muscle swelling, pulled muscles, torn muscle fibre, tendonitis, shoulder pain, tennis elbow or meniscal injuries. While **Pain relief I** provides continuous stimulation (frequency: 100 Hz / 4 Hz / 100 Hz each 10 min., pulse duration: 150 μ s at 100 Hz, 200 μ s at 4 Hz, program duration: 30 min.), **Pain relief II** offers pleasant varying stimulation (8-100 Hz, 250 μ s, 30 min.). Simply place the electrodes on or around the area of pain.

Configurable programs

U1 and **U2** are muscle build-up programs that can be individually adjusted and that include different fixed parameters, and others that can be configured according to your requirements (e.g. pulses per second, pause time, program duration). Therefore please read the chapter Programming configurable programs, page 22.

Device description

Contra-indications / application limitations

Medical circumstances under which the SaneoSPORT may not be used:

- During pregnancy, the device may never be used in the abdominal area. Current used for medical purposes can lead to miscarriage!
- When operating machinery (e.g. driving a car) or while sleeping.
- In case of metallic implants in the area of treatment.
- In the event of skin defects or defects of the mucous membrane (injuries) and acute infections of the mucous membrane, skin or subcutaneous skin.
- In case of cardiac arrhythmia or heart muscle disease .
- Where psychological limitations and/or cognitive differences are present .
- Never use in the heart region or near the carotid artery.
- Caution! The application on the thorax may lead to ventricular vibrillation.
- The device is not suitable for persons with a pacemaker and/or active implants (defibrillators)! 
- Current used for medical purposes promotes blood circulation. Treatment must be stopped and/or performed under medical supervision if:
 - A drastic change is experienced in the sensation of the current, even though the stimulation intensity was not changed.
 - Previously unknown infections or bleeding is noticed in the area of treatment prior to commencing treatment.
 - Increasing or strong skin changes are noticed, in which case a doctor must be consulted (possible allergy).

Device description

en

Scope of delivery


Item no.	Designation
ti2007	1 x SaneoSPORT Muscle Training device 1 x Instructions for use
ti6060	1 x Block battery (9 V - alkali)
ti2013	1 x body surface electrodes 50 mm x 50 mm (4 pieces)
ti2015	1 x body surface electrodes 50 mm x 90 mm (4 pieces)
ti2033	2 x Cable for connections to body surface electrodes (2-pole), 2 mm

Accessories

Item no.	Designation
	body surface electrodes:
ti2011	32 mm round (4 pieces), PZN 3888368
ti2013	50 x 50 mm (4 pieces), PZN 3888546
ti2015	50 x 90 mm (4 pieces), PZN 3888351

In Germany, selected items are also available from a pharmacy. In this case, please provide the PZN (pharmaceutical central number).

Accessory diagrams

Besides the items listed above, you can also find additional SaneoLINE products at our website:


<http://www.saneo-line.com>

Device description

Device diagrams

Overall view

en


- 1** Electrode connector socket, Channel A
- 2** Electrode connector socket, Channel B
- 3** ▲: Increase intensity
- 4** ▼: Decrease intensity
- 5** **Menu** button (select program)
- 6** **On/Off** button

Device description

Displays

All displays are listed in the following chapter Possible displays, page 18.

en

- | | | | |
|----------|---|-----------|---------------------------------|
| 1 | Battery charge status | 7 | Stimulation intensity Channel B |
| 2 | Frequency modulation active | 8 | Stimulation intensity Channel A |
| 3 | Current stimulation phase (flashing) | 9 | Invalid entry by active locking |
| 4 | Remaining treatment time in minutes:seconds | 10 | Set program |
| 5 | Locking active | 11 | Amplitude modulation active |
| 6 | Therapy session has been interrupted. | | |

Device description

Possible displays

Pxx

Muscle training program (active program)

Locking. The keypad is locked five seconds after the last entry to prevent unintentional entries. Unlocking only by means of keystroke ▼.

Invalid entry / invalid keystroke

Locking remains active.

Training session interrupted due to insufficient patient contact.

See chapter Error and problem handling (page 26), item 2.

Training time per session in minutes (**min**) and seconds.

mA

Stimulation intensity of 0 mA to 60 mA.

Burst

Mode for preventing acquired tolerance tendencies.

FM

Frequency modulation

AM

Amplitude modulation

Units of displayed values:

Hz

Frequency (pulses per second)

µs

Pulse duration (in microseconds)

sec

Stimulation and/or resting phase (in seconds)

min

Program duration (in minutes)

Battery charge status:

The battery is charged; the device is ready for use.

The battery and the device are ready for use.

Flashing: Battery must be replaced in the short term; battery charged state critical.

The current stimulation and/or resting phase is indicated by flashing. Image from left:

Pause | ▲Increase | ■Plateau | ▼Descent | Pause


Operating the SaneoSPORT

en

Preparing for application

The device is switched off; no display is present.

- ➔ Please connect the electrode cable with connector socket Channel A at the front of the device.
- ➔ Should you intend using both channels, connect an additional electrode cable with connector socket Channel B.

Preparing skin areas to be treated

Clean the skin areas to be treated of all perspiration, moisture and dirt. Ensure that the skin is dry (do not use any ointments or creams) and intact.

- ➔ Please remove the body surface electrodes from the protective packaging.
- ➔ Connect the electrodes with the electrode cable (scope of delivery).

Positioning the electrodes

Electrode application can be used to treat various indications. An electrode with two adhesive surfaces is connected to Channel A or to Channel B.

The electrode application and the program selected must match; see chapter Program summary, page 8.

- ➔ Place (stick) the electrodes as appropriate to match your program and according to the electrode application illustrated **in the cover of these instructions for use**. Observe the safety notes (page 5)!
- ➔ After treatment, the electrodes can be stuck back on the carrier sheet and stored in the protective packaging.

Operating the SaneoSPORT

Switching the SaneoSPORT on and off

Use the **On/Off** button to switch the device on and off. The software version is briefly displayed after switch-on.

en

Setting the program

Set your required program as follows (see chapter Program summary, page 8):

- ➔ Please switch the muscle training device SaneoSPORT on by means of the **On/Off** button. The *application menu* of the preset program then briefly appears in the display.
- ➔ Hold the **Menu** button down.
After approx. three seconds the *program selection menu* is displayed with the current program.
- ➔ Press the **▲** button several times until the required program is displayed.
- ➔ Briefly press the **Menu** button.
The *program selection menu* closes. The selected program is immediately displayed in the *application menu*.

Now please select the individual stimulation intensity.

Operating the SaneoSPORT

Selecting the stimulation intensity

Take note of your own perceptions!

You should clearly feel the current or stimulation intensity, but it should not be an unpleasant sensation and certainly not painful.

en

The stimulation intensity for Channel A and/or Channel B is displayed on the left and/or right-hand side of the *application menu* (unit: **mA**).

Please select the stimulation intensity for **Channel A** using the ▲/▼ buttons to the **left** of the display and for **Channel B** using the ▲/▼ buttons to the **right** of the display.

- Pressing the ▲ button briefly increases the stimulation intensity by 1 mA for the relevant channel.
- Pressing the ▼ button briefly decreases the stimulation intensity by 1 mA for the relevant channel.

Training is started.

Where the appropriate stimulation intensity was set and not changed for a period of 15 seconds, the **Lock** symbol  appears in the display. The keypad is locked.

- Where required, unlock the keypad by means of the ▼ button, e.g. to adjust the intensity.

Operating the SaneoSPORT

Programming configurable programs

U1 and **U2** are muscle build-up programs and can be changed as required.

- U1:** Basic program – amplitude modulation
(frequency and program duration configurable)
Basic program – amplitude modulation
- U2:** (rise time, plateau time, descent time, pause time
and program duration configurable)

Parameter overview

Depending on the basic program, only specific parameters can be set. Please refer to the table below for significance of the individual parameters.

The value range of the parameters is indicated (in brackets).

Parameter	Description (value range)
Hz	Training frequency in Hz (30 to 70)
	Rise time in sec (1 to 5)
	Plateau time in sec (2 to 10)
	Descent time in sec (1 to 5)
	Pause time in sec (1 to 20)
	Training time in min (5 to 40)

Operating the SaneoSPORT

en

Entering parameters

You can enter settings for training programs **U1** or **U2**. Your program configuration will then be called up as for a default program.

- ➔ Hold the **Menu** button down.
After approx. three seconds the *program selection menu* is displayed with the current program. 
- ➔ Press the **▲** button several times until the muscle build-up program to be configured is displayed.
- ➔ Press the **Menu** button to navigate to the first user-defined parameter.
The parameter value flashes and can be reduced by actuating the **▼** button or increased by actuating the **▲** button.
- ➔ Please press the **Menu** button to confirm the value selected.
The next parameter is displayed.
- ➔ Please continue as described until you reach the last parameter (**training time**).
Once you have also confirmed the training time using the **Menu** button, the *program selection menu* closes.

The *application menu* for the configured program appears in the display. Treatment with the preset parameter values can now begin.

Starting the treatment

To start the treatment, please set the stimulation intensity as described in the chapter *Selecting the stimulation intensity*, page 21.

Ending training & cleaning the electrodes

The device automatically switches off when the training time has ended. Training can be aborted at any time by using the **On/Off** button.

- ➔ Please do not remove the body surface electrodes before the device has been switched off.
- ➔


Operating the SaneoSPORT

- ➔ Use a damp, lint-free cloth to clean the electrodes after use. This will increase the durability of the electrodes.
- ➔ Stick the electrodes on the carrier sheet and store the carrier sheet in the protective packaging.

If the adhesiveness of the electrodes decreases, please use water to wet the adhesive surfaces. If this is no longer effective, please replace the electrodes.

Exchanging the battery

- ➔ Please switch the device off before exchanging the battery.

- 1** Open the battery compartment.
 - 2** Place the battery in the compartment in the correct position.
 - 3** Close the battery compartment.
- ➔ Please take note of the national statutory regulations applicable in your country in terms of environmentally friendly disposal of batteries, accessories and packaging.

Technical data and support

en

Technical data

Name and model	SaneoSPORT Muscle Training
Type	2-channel electro-stimulation system
Dimensions	105 mm x 65 mm x 26 mm
Weight	Including battery approx. 130 g
Voltage supply	Alkaline block battery (9 V) 6LR61
Classification	BF 
Current pulse shape	Biphasic rectangular pulse, charge-compensated
Relative humidity	10% to 90%
Operating temperature	0°C to 40°C
Storage temperature	0°C to 50°C
Channels	2
Current intensity	Constant Current, max. 60 mA at 500 Ω (\pm max. 5 mA)
Maximum output voltage	100 V \pm 15%
Training time	Can be adjusted from 5 min. to 40 min.
Pulse duration	150 μ s to 800 μ s
Frequency range	1 Hz to 100 Hz

Symbols

Floating application component, **Body Floating** type. The product provides protection against electrical shock in the case of permissible discharge current.

Proof of conformity in terms of Directive 93/42/EEC of the Board dated 14 June 1993 regarding medical devices. This is a medical device of class IIa (MDD).

Year of manufacture

Technical data and support

This product carries the recycling symbol in accordance with EU Directive 2002/96/EC. Disposal via domestic or residual waste is not permitted. The dealer / manufacturer is the contact person responsible for disposal.

SN

Device serial number

LOT

Lot designation

Usable until

Error and problem handling

In the event that the recommended solutions do not solve the problem, please contact Customer Services or your dealer.

1. The device does not switch on when the On/Off button is actuated.

- Ensure that the battery has been inserted correctly and contact has been established.
- Replace battery, where necessary.

2. No stimulation can be felt.

- Press the **On/Off** button to start the device.
- Check the position of the electrodes. The electrodes may not overlap.
- Gradually increase the stimulation intensity.
- Check the cable connections between the device and the electrodes.
- Please restore the adhesiveness of the electrodes by moistening.

If the adhesive surfaces of the electrodes are depleted, current distribution is no longer guaranteed.

Please replace the electrodes.