Pulse Oximeter

User Manual

Model: A320

PLEASE NOTE:

THIS MEDICAL INSTRUMENT MUST BE USED
ACCORDING TO INSTRUCTIONS TO ENSURE
ACCURATE READINGS

English	Page 1-23	Pulse Oximeter Operator's Manual
Deutsch	Page 24-46	Pulsoximeter BEDIENUNGSANLEITUNG
Français	Page 47-69	Oxymètre de pouls MANUEL DE L'UTILISATEUR
Español	Page 70-92	Oxímetro de pulso MANUAL DEL USUARIO
Italiano	Page 93-115	Pulsossimetro MANUALE D'USO

Operation quick guide

This device is intended for noninvasive spot-checking your blood oxygen saturation levels and pulse rate.

1. Install two AAA batteries into the battery compartment according to the indicated polarity + /-.

2. Insert one of your finger into the finger chamber. Finger should also be inserted completely into the chamber.

3. Turn the device on by pressing the Power button.

4. Read the Oxygen Saturation Level and Pulse Rate from the OLED Screen after a few seconds.

5. Turn the device OFF by long pressing the power button.

The device will automatically shut down for longer than 16 seconds after finger is removed from the chamber.

* Please refer to the manual for detailed operation of the Display Switching & Parameter Setting.

Table 2: Description of the number in the product diagram above

No.	Description	No.	Description
1	Power button	2	Display
3	SpO2 data	4	Battery Indicator
5	Pulse Rate data	6	Plethysmogram
7	Battery cover	8	Power off button

1 Safety

1.1 Instructions for the Safe Operation and Use of the A320 Pulse Oximeter

- Do not attempt to service the A320 Pulse Oximeter. Only qualified service personnel should attempt any needed internal servicing.
- Prolonged use or the patient's condition may require changing the sensor placement periodically. Change the sensor placement and check sensor is in contact with exposed skin, circulatory status, and correct alignment at least every 2 hours.
- Oxygen saturation (SpO2) measurements may be adversely affected in the presence of high ambient light. The sensor area should be shielded with a surgical towel, plaster or bandage if necessary.
- The following factors may cause interference to the testing accuracy of the A320 Pulse Oximeter:

- High-frequency electrosurgical equipment.
- Placement of the sensor on an extremity with a blood pressure impacted by arterial catheter, or intravascular line.
- If the patient has severe hypotension, vasoconstriction, severe anemia or hypothermia.
- If the patient is in cardiac arrest or in shock.
- Fingernail polish or false fingernails may cause inaccurate SpO2 readings.
- The device should be given at least 10 minutes to come to normal working temperature if stored in a hot or cold environment.
- The device is non-sterile and not intended to be sterilized.

1.2 Warning

The MEDICAL ELECTRONIC EQUIPMENT is suitable for home healthcare environments:

- Although the medical electronic equipment conforms to the intent of the standard EN 60601-1-2 in relation to electromagnetic compatibility, electrical equipment may produce interference. If interference is suspected, move the equipment away from the sensitive device.
- Portable and mobile RF communication equipment can affect this instrument's normal operation.
- Explosion hazard - Do not use the Pulse Oximeter in a flammable atmosphere where concentrations of flammable anesthetics or other materials may occur.

- Do not throw batteries in fire as this may cause them to explode.
- Do not attempt to recharge normal dry-cell batteries, they may leak and may cause a fire or even explode.
- Do not use the Pulse Oximeter in an MRI or CT environment.
- Do not modify this equipment without the authorization of the manufacturer.
- If this equipment is modified, appropriate inspection and testing must be conducted to ensure continued safe use of equipment.
- Do not use near active high frequency surgical equipment and the RF shielded room of an ME system for magnetic resonance imaging, where the intensity of electromagnetic disturbances is high.
- Use of this equipment adjacent to or stacked with other equipment should be avoided, because it could result in improper operation. If such use is necessary, this equipment and the other equipment should be observed to verify that they are operating normally.
- Use of accessories, transducers, and cables other than those specified or provided by the manufacturer of this equipment could result in increased electromagnetic emissions or decreased electromagnetic immunity of this equipment and result in improper operation.
- Portable RF communications equipment (including peripherals such as antenna cables and external antennas)

- should be used no closer than 30 cm (12 inches) to any part of the A320 Pulse Oximeter, including cables specified by the manufacturer. Otherwise, it may result in a degradation of the performance of this equipment.
- High-pressure sterilization cannot be used on the device.
- IF ANY: a list of all cables and maximum lengths of cables (if applicable), transducers and other accessories that are replaceable by the responsible organization and that are likely to affect compliance of the me equipment or me system with the requirements of clause 7 (emissions) and clause 8 (immunity). accessories may be specified either generically (e.g., shielded cable, load impedance) or specifically (e.g., by manufacturer and equipment or type reference).
- IF ANY: the performance of the me equipment or me system that was determined to be essential performance and a description of what the operator can expect if the essential performance is lost or degraded due to em disturbances (the defined term "essential performance" need not be used).

Cautions:

- Keep the operating environment free of dust, vibrations, corrosive, or flammable materials, and extremes of temperature and humidity.
- Do not operate the unit, if it is damp or wet because of condensation or spills. Avoid using the equipment immediately after moving it from a cold environment to a warm, humid

location.

- Never use sharp or pointed objects to operate the front-panel switches.
- The batteries must be taken out from the battery compartment, if the device will not be used for a long time.
- The device shall only be used with the battery cover closed.
- The batteries must be disposed of properly, according to local regulation after their use.
- Keep the device away from children and pets to avoid swallowing.

1.3 Definitions and Symbols

Table 3: Description of symbol

Symbol	Description	Symbol	Description
	Type BF Equipment		Batch code [1]
	Indicates the item is a medical device		Indicates a carry that contains unique device identifier information
	Information of manufacturer, including name and address		Date of manufacture [2]
			Serial No. [3]
	Temperature limitation		Informations du représentant autorisé de l'UE
	For its disposal, this product must be sent to separate collection facilities for recovery and recycling		

Symbol	Description	Symbol	Description
	Follow the instructions for use	IP22	Anti-dust & Anti-water class
<i>Caution:</i>	The information you should know to protect the equipment from possible damage.		
<i>Note:</i>	The important information you should know.		
<i>Warning</i>	The information you should know to protect patients and medical staff from possible injury.		
<i>ME</i>	Medical Electronics Equipment.		

Note: [1][2][3] Batch code, Date of manufacture and Serial No. are printed on the label on the battery cover.

2 Introduction

2.1 Intended use

The A320 Pulse Oximeter is a non-invasive device intended for spot checking of functional oxygen saturation of arterial hemoglobin (SpO₂) and pulse rate (PR). This portable device is suitable for non-professional adult patients of any color in clinical institutions and home environments.

2.2 Brief Device Description

The A320 Pulse Oximeter, is based on digital technology. The device is intended for non-invasive spot-check measurement of functional oxygen saturation of arterial hemoglobin (SpO₂). Advanced DSP algorithm [4] can minimize the influence of motion artifact and improve measurement accuracy of low perfusion [5].

2.3 Product Features

- Lightweight for easy carrying and use.
- Manually adjust the direction of interface.
- Color OLED display, simultaneous display for testing value and plethysmography [6].
- Visual & Sound reminder function. Real-time spot-checks.
- Low Battery voltage indicator.
- Automatic switch-off.
- Includes two standard AAA 1.5V alkaline batteries.

Note:[6] Plethysmograph: is an instrument for measuring changes in volume within an organ or whole body (usually resulting from fluctuations in the amount of blood or air it contains).

CAUTION: The device must not be used on children below 3 years of age, as it cannot be guaranteed to provide accurate results.

CAUTION: The Pulse Oximeter is intended only as an adjunct in patient assessment. It must be used in conjunction with other methods of assessing clinical signs and symptoms.

CAUTION: A function tester cannot be used to assess the accuracy of the Pulse Oximeter monitor or sensor.

Clinical testing is used to establish the SpO₂ accuracy. The measured arterial SpO₂ value (SpO₂) of the sensor is compared to arterial hemoglobin oxygen (SaO₂) value, determined from blood samples with a laboratory CO-oximeter. The accuracy of the sensors in comparison to the CO-oximeter samples measured over the SpO₂ range of 70-100%. Accuracy data is calculated using the root-mean-square (Arms value) for all subjects. Only about two-thirds of PULSE OXIMETER EQUIPMENT measurements can be expected to fall within $\pm$ Arms of the value measured by a CO-oximeter.

A pulse simulator shall be used to assess pulse rate accuracy. The measured pulse rate is compared to the preset pulse rate value in the simulator. Accuracy data is calculated using the root-mean-square (Arms value) for all subjects.

2.4 Expected Service Life

Is dependent on:

- The expected service life of the ME equipment;
- The expected service life of parts or accessories shipped with the ME equipment;
- Where the shelf life is less than the expected service life, the shelf life of parts or accessories shipped with the ME equipment.

3 Installation, Setup and Operation

3.1 OLED display parameter setting

When the device is in the measuring interface, press the direction button for 1 second in order to enter the menu page (figure 3.1.1 and figure 3.1.2). There are two sub-menus.

3.1.1 Reminder Setup

Press the direction button for 1 second and enter the Reminder Setup. User can adjust the setting through moving the "" symbol to the Sound Reminder, Beep, Restore or Brightness.

- Sound reminder

Press the direction button for 1 second, move the "" symbol to the "Sound Reminder", long press the direction button to turn the alarm sound on/off. If the measured value exceeds the maximum or minimum value of SPO2 or PR, there will be an alarm sound, if this function is turned on.

- Beep

Press the direction button for 1 second, move the "" symbol to the back of "Beep", long press the direction button to turn the beep sound on/off. When the beep is turned on, the sound emitted during the test indicates the pulse rate.

- Restore

With the "*" symbol next to "Restore", long press the direction button. This will change to "OK", which restores the device to factory settings.

- Demo

With the "*" symbol next to "Demo", long press the direction button to turn Demo mode on/off.

- Brightness

With the "*" symbol next to "Brightness", long press the direction button to change the Brightness value on a scale from 1 to 5.

3.1.2 Limit Value Setting

With the "*" symbol next to "Reminder Setup", long press the direction button until the "Limit Setup" menu appears (figure 3.3.2). You can then press the direction button to select the items. Press the direction button for 1 second to change to the data you need.

Remind Setup	*
Sound Reminder	on
Beep	off
Demo	on
Restore	OK
Brightness	4
Exit	

Figure 3.3.1

Limit Setup	*
SpO2 Hi	100
SpO2 Lo	94
PR Hi	130
PR Lo	50
+/-	+
Exit	

Figure 3.3.2

On the "Limit Setup" menu page (figure 3.3.2), with the "*" symbol next to "+/-", press the direction button for 1 second to change the "+" to "-" or change the "-" to "+".

When "+" shows on the right side, by pressing the direction button for 1 second and moving the "*" next to the Spo2 Hi or PR Hi setting, you can increase the value to a higher value (until it reaches the

maximum).

When “-” shows on the right side, by press the direction button for 1 second and moving the “*” next to the Spo2 Lo or PR Lo value setting, you can reduce the value to a lower value (until it reaches the minimum).

3.2 Operation

3.2.1 Installing batteries


Place two AAA batteries into the battery compartment taking care to align the polarities, then close the cover.

WARNING: Do not attempt to recharge normal alkaline batteries, as they may leak and may cause a fire or even explode.

3.2.2 Turning on and applying the Pulse Oximeter

Put one of your fingers into the rubber opening of the A320 Pulse Oximeter with nail facing upward, then release the clamp. Press the power button to turn the A320 Pulse Oximeter on. The oximeter will automatically turn off, if there is no finger in the device for more than 16 ± 2 seconds.

3.2.3 Read data from display screen

Table 4: Display Description of OLED

Display	OLED Screen		
Model	A320		
Description	The screen display can scroll through four directions with six different display modes by pressing the power button.		
Show	%SpO₂ PR bpm 98 65 PI% 5.5 %SpO₂ PR bpm 98 65 PI% 5.5 %SpO₂ PR bpm 98 65 PI% 5.5 %SpO₂ PR bpm 65 86 PI% 5.5 %SpO₂ PR bpm 98 65 PI% 5.5 %SpO₂ PR bpm 98 65 PI% 5.5 		
Note	1. When the battery power is at its lowest level, the battery symbol will be shown, reminding users to replace the batteries. 2. The plethymogram can be considered accurate if the wave symbol is fluctuating regularly.		

4 Cleaning and Disinfection

4.1 Cleaning

Switch off the power and take out the batteries before cleaning. Keep the exterior surface of the device clean and free of dust and dirt. Clean the exterior surface (display screen included) of the unit with a soft dry cloth. Use 75% medical alcohol to clean the surface by applying a small amount with a dry cloth to avoid the alcohol entering the device.

4.2 Disinfection

Disinfect the machine after each use, if multiple patients use the machine, e.g., in a hospital.

Use 75% medical alcohol to clean the surface that was in contact with the patient.

CAUTION: Do not use strong solvent, e.g., acetone.

CAUTION: Never use an abrasive such as steel wool or metal polish.

CAUTION: Do not allow any liquid to enter the product and do not immerse any part of the device in liquid.

CAUTION: Avoid pouring liquid on the device while cleaning.

CAUTION: Do not leave cleaning solution on the surface of the device.

5 Troubleshooting and Maintenance

5.1 Maintenance

Replace the batteries in a timely manner, if the battery indication is low. Clean the surface of the oximeter before it is used for diagnosis of patients. Remove the batteries from the battery compartment if the oximeter will not be operated for a long time. It is best to store the product where the ambient temperature is between -25°C and 55°C and humidity is 15% to 93%. Regular inspection is recommended in order to make sure that no obvious damage is present that may affect the safety and performance of the device. Do not expose the device to flammable substances, high or low temperatures or humidity levels outside those in the operation conditions.

5.2 Troubleshooting

Table 5: Troubleshooting

Problems	Possible Reason	Resolutions
Oxyhemoglobin or heart rate cannot be shown.	1. Finger is not inserted correctly. 2. Patient's perfusion is too low to be measured.	1. Retry by inserting the finger thoroughly. 2. Try a few more times to make sure there is no problem with the product itself. Otherwise seek medical help for exact diagnosis.
Oxyhemoglobin or heart rate is unstable.	1. Finger might not be inserted correctly. 2. Finger is trembling, or patient's body is moving.	1. Retry by inserting the finger thoroughly. 2. Try to help the patient keep calm and still.

Oxyhemoglobin or heart rate is outside of the standard range.	1. Finger might not be inserted correctly. 2. Patient's SPO2 & PR is abnormal.	1. Retry by inserting the finger thoroughly. 2. Seek medical help for further examination.
The Oximeter cannot be turned on.	1. Batteries may need replacement. 2. Batteries might be installed incorrectly. 3. The oximeter might be damaged.	1. Please replace batteries. 2. Please reinstall the batteries. 3. Please contact your point of purchase.
The screen suddenly turns off.	1. The device automatically turns off, if there is no signal detected for more than 16 seconds. 2. Batteries may need replacement.	1. This is a normal process. 2. Replace the batteries.

6 Specification

Table 6: Specification of A320 pulse oximeter

Device Name	Pulse Oximeter
Dimensions (L×W×H)	73*37*38 mm
Weight	Approx.50g - 60g (including 2 × AAA battery)
Anti-electric Shock Type	Internally powered equipment
Anti-electric Shock Equipment Degree	Type BF
EMC Type	Group I Class B
Enclosure Degree of ingress protection	IP22
Internal Power:	2xAAA 1.5v alkaline battery
Power Consumption	Below 45mA
Screen	0.96' OLED
SpO2 Display	35-100%
Pulse rate Display	30-250 BPM
Resolution	SpO2: 1%
	Pulse rate: 1BPM
Measure Accuracy	SpO2: ±3% (70%-100%) Unspecified (<70%)
	PR: ±2BPM

Data averaging and other signal processing	8s
Data Update Period	1s
Operating Environment	Temperature: +5°C to +40°C Humidity: 15% to 93% non-condensing Air Pressure: 70Kpa-106Kpa
Storage & Transport Environment	Temperature: -25°C to +55°C Humidity: 15% to 93% non-condensing Air Pressure: 70Kpa-106Kpa

7 Clinical SpO2 Accuracy

The below table shows statistic distribution of an invasive controlled desaturation study, which guided by ISO80601-2-61, Annex EE, Guideline for evaluating and documenting SpO2 Accuracy in human subjects.

The statistic distribution displayed the accuracy distribution between the range of 70% ~ 100%, which may helpful to user.

Table 7: Clinical SpO2 Accuracy

Device	Item	100-70%	100 - 90%	90- 80%	80 - 70%
A320	Bs	1.0	0.80	1.26	0.86
	Sres	1.94	0.86	1.61	2.91
	Arms	2.18	1.17	2	2.99

The below is the Bland-Altman graphical plot of samples from invasive controlled desaturation study.


Figure 7.1-Bland & Altman


Figure 7.2-Linear regression for SpO2 vs SaO2

8 Disposal

Consider the applicable regulations when disposing of the A320 Pulse Oximeter and batteries. The pulse oximeter must not be disposed of in the domestic waste. All users are obliged to hand in all electrical or electronic devices, regardless of whether they contain toxic substances, at a municipal or commercial collection point so that they can be disposed of in an environmentally acceptable manner. Please remove the batteries before disposing of the pulse oximeter. Do not dispose of old batteries with your household waste, but at a battery collection station at a recycling site or in a shop.

9 Certificate of guarantee

Aeon guarantees the A320 Pulse Oximeter against any manufacturing defect for one year from the date of purchase, if it is returned to the dealer from which it was purchased in original packaging and with the purchase receipt.

During this period the unit will be repaired or replaced free of charge, if the fault is due to defective design or assembly.

This guarantee does not cover any damage or defect caused by improper handling resulting from use that is not in compliance with these instructions or from unauthorized attempts to repair it.

10 Manufacturer's EMC Declaration

- 1) All necessary instructions for maintaining BASIC SAFETY and ESSENTIAL PERFORMANCE with regard to electromagnetic disturbances for the expected service life.
- 2) Guidance and manufacturer's declaration - Electromagnetic Emissions and Immunity.

Table 8 - Guidance and manufacturer's declaration - electromagnetic emissions	
Emissions test	Compliance
RF emissions CISPR 11	Group 1
RF emissions CISPR 11	Class B
Harmonic emissions IEC 61000-3-2	Not application
Voltage fluctuations/ flicker emissions IEC 61000-3-3	Not application

Table 9 - Guidance and manufacturer's declaration - Electromagnetic Immunity		
Immunity Test	IEC 60601-1-2 Test level	Compliance level
Electrostatic discharge (ESD) IEC 61000-4-2	± 8 kV contact ± 2 kV, ± 4 kV, ± 8 kV, ± 15 kV air	± 8 kV contact ± 2 kV, ± 4 kV, ± 8 kV, ± 15 kV air
Electrical fast transient/burst IEC 61000-4-4	Power supply lines: ± 2 kV input/output lines: ± 1 kV	Not application
Surge IEC 61000-4-5	line(s) to line(s): ± 1 kV. line(s) to earth: ± 2 kV. 100 kHz repetition frequency	Not application

Voltage dips, short interruptions, and voltage variations on power supply input lines IEC 61000-4-11	0% 0.5 cycle At 0°, 45°, 90°, 135°, 180°, 225°, 270° and 315° 0% 1 cycle And 70% 25/30 cycles Single phase: at 0 0% 300 cycle	Not application
Power frequency magnetic field IEC 61000-4-8	30 A/m 50Hz/60Hz	30 A/m 50Hz/60Hz
Conducted RF IEC61000-4-6	150KHz to 80MHz: 3Vrms 6Vrms (in ISM and amateur radio bands) 80% Am at 1kHz	Not application
Radiated RF IEC61000-4-3	10 V/m 80 MHz – 2,7 GHz 80 % AM at 1 kHz	10 V/m 80 MHz -2,7 GHz 80 % AM at 1 kHz
NOTE UT is the a.c. mains voltage prior to application of the test level.		

Table 10 - Guidance and manufacturer's declaration - Electromagnetic Immunity

Radiated RF IEC61000-4-3 (Test specifications for	Test Frequency (MHz)	Band (MHz)	Service	Modulation	Modulation (W)	Distance (m)	Immunity test Level (V/m)
	385	380-390	TETRA 400	Pulse Modulation 18 Hz	1,8	0.3	27
	450	430-470	GMRS 460, FRS 460	FM ± 5 kHz deviation 1 kHz sine	2	0.3	28
	710 745 780	704-787	LTE Band 13, 17	Pulse Modulation 217 Hz	0,2	0.3	9

ENCLOSURE PORT IMMUNITY to RF wireless communications equipment)	810 870 930	800- 960	GSM 800/900, TETRA 800, iDEN 820, CDMA 850, LTE Band 5	Pulse Modulation 18 Hz	2	0.3	28
	1720 1845 1970	1700- 1990	GSM 1800; CDMA 1900; GSM 1900; DECT; LTE Band 1,3,4, 25; UMTS	Pulse Modulation 217 Hz	2	0.3	28
	2450	2400- 2570	Bluetooth, WLAN, 802.11 b/g/n, RFID 2450, LTE Band 7	Pulse Modulation 217 Hz	2	0.3	28
	5240 5500 5785	5100- 5800	WLAN 802.11 a/n	Pulse Modulation 217 Hz	0,2	0.3	9

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